

Signal™ Instrument Transport Bags/Labels.

Protect staff, patients, and visitors from biohazard risks.

Place a high-level disinfected instrument in the Signal Instrument Transport Bag and seal the bag using the low-tack re-sealable closure. Use the included HLD label (order #ITB1218-H) to record the date the package was prepared, and the preparer's name. Once the procedure is over, and the instrument is to be returned to decontamination, remove the HLD label to expose the biohazard legend. Place the soiled instrument back in the bag and re-seal the closure.



Close-up of the HLD Label

Instrument Transport Bag with Removable HLD Label

- Labels with low-tack adhesive are easy to remove even while wearing bulky gloves.
- Clear easy-to-read color-coded green high-level-disinfected labels.
- The low-tack adhesive won't leave residue behind.
- The color-coded red Biohazard Legend is pre-printed on the bag.

Bag with HLD labels included:

ITB1218-H

Signal Instrument Transport Bag, resealing, with biohazard legend - 200 HLD Labels included.

Size: 12" x 18"

Units/Case: 200 Bags/Labels

Bag without HLD labels:

ITB1218

Signal Instrument Transport Bag, resealing, with biohazard legend.

Size: 12" x 18"

Units/Case: 200 Bags



Cygnus Medical 800.990.7489

965 West Main Street, Branford, CT 06405

Email: sales@cygnusmedical.com | www.cygnusmedical.com

SIGNAL 
A Label Revolution