

Hoping instrument jaws stay open during sterilization cycles may not be enough.

Achieve  of mind with the Gator Instrument Clip.

Exposing instrument surface area during sterilization is critical in order to ensure sterility. Gator Instrument Clips from Cygnus Medical help to keep instrument jaws open during the sterilization cycle.

Peel pouch instruments must be sterilized on their side edges. In this position gravity pulls the instruments closed. Gator Instrument Clips prevent this from happening.



- Latex-Free
- Single-Use
- Multiple Sizes
- Radiopaque



Cygnus Medical

965 West Main Street, Branford, CT 06405
800.990.7489 www.cygnusmedical.com



PATENT PENDING