

SIGNAL™

A Label Revolution

Instrument Transport Bag with Removable HLD Label

Clear communication protects staff, patients, and visitors from biohazard risks.



Place a High-Level Disinfected Label over the biohazard legend printed on the Transport Bag to identify high-level disinfected contents. Remove the HLD label when contents become soiled to uncover the biohazard legend.

- Easy to apply and remove even while wearing gloves.
- The low-tack adhesive won't leave residue behind.
- Compostable (ASTM D6400 standards) and made of a biobased material.

Instrument Transport Bag

Single-Use & Latex-Free | Order#: ITB2020



- This film is 100% certified compostable for both industrial composting facilities or in a home compost pile.
- Adheres to the ASTM D6400 standards on compostability.
- Composed of biobased material.



965 West Main Street, Branford, CT 06405 | 800.990.7489

TO REMOVE LABEL - PULL DOWN HERE

SIGNAL™
A Label Revolution

High-Level Disinfected

Date: _____

Prepared by: _____

800.990.7489 | www.cygnusmedical.com
ORDER#: HLD500

ITB1218

Signal Instrument Transport Bag, compostable, resealing, with biohazard legend.

Size: 12" x 18"

Units/Case: 200 Bags

ITB1218-H

Signal Instrument Transport Bag, compostable, resealing, with biohazard legend - HLD Label pre-applied.

Size: 12" x 18"

Units/Case: 200 Bags/Labels

ITB2022

Signal Instrument Transport Bag, compostable, resealing, with biohazard legend.

Size: 20" x 22"

Units/Case: 200 Bags

ITB2022-H

Signal Instrument Transport Bag, compostable, resealing, with biohazard legend - HLD Label pre-applied.

Size: 20" x 22"

Units/Case: 200 Bags/Labels

HLD500

Signal High-Level Disinfected Label with low-tack adhesive.

Size: 4" x 4"

Units/Case: 500 Labels



Cygnus Medical 800.990.7489

965 West Main Street, Branford, CT 06405

Email: sales@cygnusmedical.com | www.cygnusmedical.com