



POINT-OF-USE

Canteen™

PRE-SATURATED

Surgical Instrument Wipers

ORDER#: CW150
Pack of Five - 7" x 8" Wipers

CYGNUS
MEDICAL

965 West Main Street
Branford, CT 06405
www.cygnusmedical.com

- Canteen Wipers are Pre-Saturated with Deionized Water.
- Removes Gross Contamination on Surgical Instruments.
- Non-Sterile, Post-Procedure, Point-of-Use Cleaning Only.
- Clean Wipers Can Be Used to Drape Over Instruments During Transport to Maintain Moisture.
- Single-Use & Latex-Free.

CLEANING DILEMMA | Your surgical instruments need to be rinsed, no Utility Sink is nearby, and the Hand Washing Sink is NOT an option?

You'll need your canteen.



- Pre-Saturated with 94% Deionized Water and 6% IPA*
- Removes gross contamination on surgical instruments.
- Non-Sterile, post-procedure, point-of-use cleaning only.
- Clean wipers can be draped over instruments during transport to maintain moisture.
- Single-Use & Latex-Free.

*IPA used as a stabilizer only.

Order#: CW150
Size: 7" x 8" Pre-Saturated Wipers
Case Quantity: 750 Wipers Total (5 Wipers to a Package).

POINT-OF-USE

Canteen™

PRE-SATURATED

Surgical Instrument Wipers



965 West Main Street, Branford, CT 06405
www.cygnusmedical.com | 800.990.7489