


# Together, but Separate.

## DuPlex™ Dual-Cavity Instrument Transport Bag




**HEAT-SEALED STRIP**  
Separates the **BIOHAZARD**  
compartment from the  
**CORD** compartment.

### A Transport Solution

Not all parts of internal probes require high-level disinfection (HLD). Cords and connectors need to be kept separate from the internal transducer that has been exposed to gross contamination.

DuPlex™ Dual-Cavity Instrument Transport Bags isolate dirty portions of instruments and prevent the cross-contamination of sections that are not high-level disinfected.



DUAL-CAVITY INSTRUMENT TRANSPORT BAG

**Order#: DCB2017**

**Size: 17" x 20"**

**Case Quantity: 100 Bags**



965 West Main Street, Branford, CT 06405  
www.cygnusmedical.com | 800.990.7489